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

50 52 T

EN

For more detailed information, refer to the meter documentation.

Caution
 ⚠ Personal injury hazard. Glass components can break. Handle with care to prevent cuts.
 ⚠ Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.

Introduction
 This probe is a combination gel-filled electrode with a built-in temperature sensor for pH measurements in low-conductivity, dirty, viscous and high temperature samples.

Technical information

Range	0 to 14 pH
Junction	PTFE
Reference	Ag/AgCl
Temperature	Pt 1000
Minimum depth	12 mm (0.5 in.)
Operating temp	0–100 °C (32–212 °F)
Connector type	MP-5

Preparation
 Install the protective shroud.
Connect the probe
 ⚠ Do not turn or twist the probe connector.

1 Calibration
 ⚠ Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.

For 1 or 2-point pH calibration: push ✓ after 1 or 2 standards.
2 Sample measurement
 Continuous: push ✓ 2 times.
 Change the parameter: push and hold ✓ during stabilization.

Maintenance and storage
 Cleaning: soak the probe in a cleaning solution. Rinse with deionized water.
 Hydration: soak the dry probe in a dilute HCl solution for several hours.
 Storage: put the probe in the protector with KCl 3 M solution.

Troubleshooting
 Be sure that the correct probe type is used. Do not hold the tubes by hand to prevent temperature changes.

Warranty
 6 months for manufacturing faults. This warranty does not cover improper use or wear.

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