ISE Na⁺ 9650

User Manual

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


www.hach.com

EN

ISE Na⁺ 9650
Refer to the meter manual for more information.

Caution
 ⚠ *Personal injury hazard. Glass components can break. Handle with care to prevent cuts.*
 ⚠ *Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.*

Introduction
This probe measures sodium concentration in aqueous solutions. Use with a reference electrode.

Technical information

Range	0.02 mg/L to 23 g/L Na ⁺
pH range	8 to 11
Linear range	1 mg/L to 23 g/L Na ⁺
Slope	-59 mV/pNa, 90 to 110% at 25 °C (77 °F)
Ionic strength adjuster (ISA)	Hach ISA 4451569
Operating temp	0 to 60 °C (32 to 140 °F)
Interferences	K ⁺ (5), NH ₄ ⁺ (50), H ⁺ (0.001), Li ⁺ (1)

Preparation
Soak the electrode in a 0.01 M NaCl solution for 20 to 30 minutes.

1 Connect the probe
2 Calibration
Add 1 ISA pillow per 25 mL of standard.
3 Measurement
Add 1 ISA pillow per 25 mL of sample.

Maintenance and storage
Cleaning: rinse with distilled water between each standard and each sample.
Storage: clean and dry the probe for storage.

Troubleshooting
The ionic strength and temperature of samples and standards must be the same. Make sure the stir speed is the same for all standards and samples.
Damage or deposits on the membrane change the response of the electrode.

Warranty
6 months for manufacturing faults only. Damage from use is not covered.

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